


Moving Head type Spectroscopic Ellipsometer

- Ultra-large area mapping spectroscopic ellipsometer
 - Easy to Use
 - Ultra thin film Sensitivity
 - Best for Routine Measurement of Dielectrics



- Moving head type (patented pending)
 - > Unlimited sample size
 - ※ Conventional mapping ellipsometers : moving sample type (limited)
- Fast spectroscopic data
 - > 5 sec for full spectra of $\{\Delta, \Psi\}$ over 370~830nm (UV optional)
 - > Multi-head type is possible for in line application
- Easiest operation in the world
 - > Self alignment, self calibration
 - > maintenance-free (except lamp)
- User-friendly operation and analysis software
- Other features
 - > fixed angle of incidence : $70^\circ \pm 0.5^\circ$ (customize)
 - > computer (window 2000, XP) : optional